

SCANLAN® EASY-TAG™ Tracking System

RAPID TURNOVER

ROOM # 211 TIME NEEDED 07:30
SURGEON Anderson

SCANLAN ★

INTERNATIONAL

Surgical products designed and manufactured by the Scanlan family since 1921

SCANLAN® EASY-TAG™ Tracking System

- **Sterilizable** - EO, Steam and Gas Plasma
- **Versatile Surface** - allows placement of pre-printed labels or can be written on both sides of tag
- **No Adhesives** - quick, easy and secure
- **Easy to Apply** - the flexible handle loops around instruments and tray handles
- **Easy to Remove** - Tear-Away™ handle provides easy separation after use
- **Color Coded** - for fast and easy identification
- **Not made with natural rubber latex**
- 250 tags per spool

EASY-TAG™ Tracking System

2.5" x 4.75" (63.5 mm x 120.65 mm)

Cat. No. Description

1001-904	EASY-TAG™, Black
1001-905	EASY-TAG™, Red
1001-906	EASY-TAG™, Yellow
1001-907	EASY-TAG™, Blue
1001-908	EASY-TAG™, Green
1001-909	EASY-TAG™, Pink
1001-910	EASY-TAG™, Orange

EASY-TAG™ Tracking System Small

2.5" x 2.875" (63.5 mm x 73.03 mm)

Cat. No. Description

1001-914	EASY-TAG™, White
1001-915	EASY-TAG™, Red
1001-916	EASY-TAG™, Purple
1001-917	EASY-TAG™, Light Blue
1001-918	EASY-TAG™, Lime Green
1001-919	EASY-TAG™, Hot Pink
1001-920	EASY-TAG™, Orange

Use your Smart Phone and QR Reader software
to photograph the code below to view a
**Video Demonstration of the
EASY-TAG™ Tracking System**



EASY-TAG™ SMALL





SCANLAN® EASY-TAG™ Pre-Printed

250 tags per spool

EASY-TAG™ Pre-Printed

2.5" x 4.75" (63.5 mm x 120.65 mm)

Cat. No.	Description
1001-921	EASY-TAG™, Blue Endoscope Cleaned
1001-922	EASY-TAG™, Green Endoscope Cleaned Month / Date
1001-923	EASY-TAG™, Red Rapid Turnover
1001-924	EASY-TAG™, Red Repair
1001-925	EASY-TAG™, Green Clean Ready To Use
1001-926	EASY-TAG™, Black Loaner
1001-927	EASY-TAG™, Red Biohazard

1001-921

1001-922

1001-923

1001-924

1001-925

1001-926

1001-927



SCANLAN® EASY-TAG™ Custom Printing

Customize your EASY-TAG™ Tracking System to meet the unique needs within your hospital. Simply define your message, and we do the rest! Your custom print EASY-TAG™ will be developed and printed to your specific needs.

250 tags per spool

Call your Sales Rep at 800-328-9458 to customize your Easy-Tag™.



Custom EASY-TAG™

2.5" x 4.75" (63.5 mm x 120.65 mm)

Cat. No. Description

1001-904C	EASY-TAG™, Black
1001-905C	EASY-TAG™, Red
1001-906C	EASY-TAG™, Yellow
1001-907C	EASY-TAG™, Blue
1001-908C	EASY-TAG™, Green
1001-909C	EASY-TAG™, Pink
1001-910C	EASY-TAG™, Orange

Not made with natural rubber latex
Custom made at time of order. Firm sale.

Custom EASY-TAG™ Small

2.5" x 2.875" (63.5 mm x 73.03 mm)

Cat. No. Description

1001-914C	EASY-TAG™, White
1001-915C	EASY-TAG™, Red
1001-916C	EASY-TAG™, Purple
1001-917C	EASY-TAG™, Light Blue
1001-918C	EASY-TAG™, Lime Green
1001-919C	EASY-TAG™, Hot Pink
1001-920C	EASY-TAG™, Orange

Not made with natural rubber latex
Custom made at time of order. Firm sale.

SCANLAN® EASY-TAG™ Plus

Patented flexible/elastic handle technology of the Easy-Tag Tracking System, PLUS an adhesive back allowing the user to attach a machine printout, note, or message.

- **Easy to Apply** - the flexible handle loops over endoscope knobs, around instruments, tray handles and carts.
- **Adhesive Back** - allows placement of a machine print out, note or message to the elastic handle
- **Not made with natural rubber latex**
- 500 tags per spool



EASY-TAG™ Plus

Cat. No. Description

1001-940	Easy-Tag Plus, Silver *
-----------------	-------------------------

* Not intended for sterilization



For more information or to
place an order,
please contact your local
Scanlan representative, or call...

U.S. & Canada: 800-328-9458
International: 651-298-0997
FAX: 651-298-0018
www.scanlaninternational.com



SCANLAN ★

INTERNATIONAL

Surgical products designed and manufactured by the Scanlan family since 1921
ONE SCANLAN PLAZA • SAINT PAUL, MINNESOTA 55107 • U.S.A.
800-328-9458 • 651-298-0997 • FAX 651-298-0018 • www.scanlaninternational.com

MK-B-EZTS Rev. B
© 2015 Scanlan International, Inc.