

SURG-I-BAND®

Color Coding Data Matrix Barcode



INTERNATIONAL

Surgical products designed and manufactured by the Scanlan family since 1921

SURG-I-BAND®

Color Coding Data Matrix Barcode
















Offering precise instrument / device tracking

Surg-I-Band® Data Matrix offers hospitals an excellent system to organize and track surgical instrumentation and devices. Combining the benefits of Surg-I-Band® color coding with the technology of bar coding, Surg-I-Band® Data Matrix provides a fast, easy and precise organizational system.

Advantages include...

- **Economy** – Each spool contains 200 unique data matrix barcodes
- **14 colors** – For fast and easy identification
- **Safety** – Clinically tested for steam, gas and flash sterilization
- **Trouble Free** – Easy application and removal without damaging instrument
- **Versatility** – Can be used on most sizes and designs of instruments
- Compatible with instrument tracking software and high-density data matrix ECC 200 scanners
- May be used in combination with Surg-I-Band® Dispenser, Remover and Carousel



Description		Data Matrix	Censitrac™ Compatible
Red			1001-1084 1001-1098
Yellow			1001-1085 1001-1099
Orange			1001-1086 1001-1100
Blue			1001-1087 1001-1101
Green			1001-1088 1001-1102
Brown			1001-1089 1001-1103
White			1001-1090 1001-1104
Black			1001-1091 1001-1105
Glo-Orange			1001-1092 1001-1106
Purple			1001-1093 1001-1107
Lime Green			1001-1094 1001-1108
Grey			1001-1095 1001-1109
Pink			1001-1096 1001-1110
Light Blue			1001-1097 1001-1111



Application Video

Please visit the Surg-I-Band® Color Coding page at www.scanlaninternational.com or scan the QR Code.

For more information or to place an order,
please contact your local Scanlan representative, or call...

U.S. & Canada: 800-328-9458

International: 651-298-0997

FAX: 651-298-0018

www.scanlaninternational.com



INTERNATIONAL

Surgical products designed and manufactured by the Scanlan family since 1921

ONE SCANLAN PLAZA • SAINT PAUL, MINNESOTA 55107 • U.S.A.
800-328-9458 • 651-298-0997 • FAX 651-298-0018 • www.scanlaninternational.com

MK-B-SIBDM Rev. B
© 2014 Scanlan International, Inc.